

SLC46A2 Antibody, FITC conjugated



PACO41008

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May act as a transporter.

Gene ID:

SLC46A2

Uniprot

Q9BY10

Synonyms:

Thymic stromal cotransporter homolog (Solute carrier family 46 member 2), SLC46A2, TSCOT

Immunogen:

Recombinant Human Thymic stromal cotransporter homolog protein (45-78AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A