

C22orf23 Antibody, FITC conjugated



PACO41020

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C22orf23

Uniprot

Q9BZE7

Synonyms:

UPF0193 protein EVG1, C22orf23

Immunogen:

Recombinant Human UPF0193 protein EVG1 protein (1-217AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A