

WDR61 Antibody, Biotin conjugated



PACO41037

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

WDR61

Uniprot

Q9GZS3

Synonyms:

WD repeat-containing protein 61 (Meiotic recombination REC14 protein homolog) (SKI8 homolog) (Ski8) [Cleaved into: WD repeat-containing protein 61, N-terminally processed], WDR61

Immunogen:

Recombinant Human WD repeat-containing protein 61 protein (2-305AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A