

SDCBP2 Antibody



PACO41054

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytoplasm, extracellular exosome, plasma membrane, protein C-terminus binding, protein heterodimerization activity, protein homodimerization activity, intracellular signal transduction, intracellular transport, nervous system development.

Gene ID:

SDCBP2

Uniprot

Q9H190

Synonyms:

Syntenin-2 (Similar to TACIP18) (SITAC) (Syndecan-binding protein 2), SDCBP2, SITAC18


Immunogen:

Recombinant Human Syntenin-2 protein (11-207AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO41054 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO41054 at dilution of 1:100.