

SLC25A32 Antibody, HRP conjugated



PACO41059

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transports folate across the inner membranes of mitochondria.

Gene ID:

SLC25A32

Uniprot

Q9H2D1

Synonyms:

Mitochondrial folate transporter/carrier (Solute carrier family 25 member 32), SLC25A32, MFT MFTC

Immunogen:

Recombinant Human Mitochondrial folate transporter/carrier protein (44-88AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A