

# SLC25A32 Antibody, FITC conjugated



PACO41060

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transports folate across the inner membranes of mitochondria.

**Gene ID:**

SLC25A32

**Uniprot**

Q9H2D1

**Synonyms:**

Mitochondrial folate transporter/carrier (Solute carrier family 25 member 32), SLC25A32, MFT MFTC

**Immunogen:**

Recombinant Human Mitochondrial folate transporter/carrier protein (44-88AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A