

FBXO44 Antibody



PACO41074

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

Gene ID:

FBXO44

Uniprot

Q9H4M3

Synonyms:

F-box only protein 44 (F-box protein FBX30) (F-box/G-domain protein 3), FBXO44, FBG3, FBX30, FBX44, FBX6A, FBXO6A

Immunogen:

Recombinant Human F-box only protein 44 protein (1-224AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A