

ARMC7 Antibody



PACO41086

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytoplasm.

Gene ID:

ARMC7

Uniprot

Q9H6L4

Synonyms:

Armadillo repeat-containing protein 7, ARMC7


Immunogen:

Recombinant Human Armadillo repeat-containing protein 7 protein (1-198AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO41086 at dilution of 1:100.