

# DDX31 Antibody, HRP conjugated



PACO41099

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable ATP-dependent RNA helicase. Plays a role in ribosome biogenesis and TP53/p53 regulation through its interaction with NPM1.

**Gene ID:**

DDX31

**Uniprot**

Q9H8H2

**Synonyms:**

Probable ATP-dependent RNA helicase DDX31 (EC 3.6.4.13) (DEAD box protein 31) (Helicain), DDX31

**Immunogen:**

Recombinant Human Probable ATP-dependent RNA helicase DDX31 protein (560-851AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A