

FAM49B Antibody



PACO41174

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FAM49B

Uniprot

Q9NUQ9

Synonyms:

Protein FAM49B (L1), FAM49B

Immunogen:

Recombinant Human Protein FAM49B protein (1-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A