

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Small GTPase-like component of the intraflagellar transport (IFT) complex B that promotes the exit of the BBSome complex from cilia via its interaction with ARL6. Not involved in entry of the BBSome complex into cilium. Prevents aggregation of GTP-free ARL6. Required for hedgehog signaling. Forms a subcomplex within the IFT complex B with IFT25.

Gene ID:

IFT27

Uniprot

Q9BW83

Synonyms:

Intraflagellar transport protein 27 homolog (Putative GTP-binding protein RAY-like) (Rab-like protein 4), IFT27, RABL4 RAYL


Immunogen:

Recombinant Human Intraflagellar transport protein 27 homolog protein (1-186AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO41262 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO41262 at dilution of 1:100.