

# STRADB Antibody, HRP conjugated



PACO41287

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Pseudokinase which, in complex with CAB39/MO25 (CAB39/MO25alpha or CAB39L/MO25beta), binds to and activates STK11/LKB1. Adopts a closed conformation typical of active protein kinases and binds STK11/LKB1 as a pseudosubstrate, promoting conformational change of STK11/LKB1 in an active conformation.

**Gene ID:**

STRADB

**Uniprot**

Q9C0K7

**Synonyms:**

STE20-related kinase adapter protein beta (STRAD beta) (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 2 protein) (CAL5-21) (ILP-interacting protein) (Pseudokinase ALS2CR2), STRADB, ALS2CR2 ILPIP

**Immunogen:**

Recombinant Human STE20-related kinase adapter protein &beta; protein (212-418AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A