

MIS12 Antibody, Biotin conjugated



PACO41313

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Part of the MIS12 complex which is required for normal chromosome alignment and segregation and for kinetochore formation during mitosis.

Gene ID:

MIS12

Uniprot

Q9H081

Synonyms:

Protein MIS12 homolog, MIS12

Immunogen:

Recombinant Human Protein MIS12 homolog protein (1-305AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A