

DOCK5 Antibody



PACO41346

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Guanine nucleotide exchange factor (GEF) for Rho and Rac. GEF proteins activate small GTPases by exchanging bound GDP for free GTP.

Gene ID:

DOCK5

Uniprot

Q9H7D0

Synonyms:

Dedicator of cytokinesis protein 5, DOCK5

Immunogen:

Recombinant Human Dedicator of cytokinesis protein 5 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A