

RMND5A Antibody, HRP conjugated



PACO41351

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

RMND5A

Uniprot

Q9H871

Synonyms:

Protein RMD5 homolog A, RMND5A

Immunogen:

Recombinant Human Protein RMD5 homolog A protein (1-391AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A