

PACO41366

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:**

Neuronal cell surface protein that may be involved in cell recognition and cell adhesion. May play a role in angiogenesis.

**Gene ID:**

NRXN3

**Uniprot**

Q9HDB5

**Synonyms:**

Neurexin-3-beta (Neurexin III-beta) [Cleaved into: Neurexin-3-beta, soluble form; Neurexin-3-beta, C-terminal fragment (NRXN3-CTF)], NRXN3, KIAA0743

**Immunogen:**

Recombinant Human Neurexin-3-&beta; protein (201-432AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western blot. All lanes: NRXN3 antibody at 2 $\mu$ g/ml + Mouse heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 70, 48, 44, 51 kDa. Observed band size: 70 kDa.