

SOHLH2 Antibody



PACO41450

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Probable transcription factor, which may be involved in spermatogenesis and oogenesis.

Gene ID:

SOHLH2

Uniprot

Q9NX45

Synonyms:

Spermatogenesis- and oogenesis-specific basic helix-loop-helix-containing protein 2, SOHLH2, TEB1


Immunogen:

Recombinant Human Spermatogenesis- and oogenesis-specific basic helix-loop-helix-containing protein 2 protein (1-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO41450 at dilution of 1:100.