

Prok2 Antibody, Biotin conjugated



PACO41505

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May function as an output molecule from the suprachiasmatic nucleus (SCN) that transmits behavioral circadian rhythm. May also function locally within the SCN to synchronize output. Potently contracts gastrointestinal (GI) smooth muscle.

Gene ID:

Prok2

Uniprot

Q9QXU7

Synonyms:

Prokineticin-2 (PK2) (Protein Bv8 homolog), Prok2, Bv8

Immunogen:

Recombinant Mouse Prokineticin-2 protein (27-128AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A