

PISD Antibody



PACO41514

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PISD

Uniprot

Q9UG56

Synonyms:

Phosphatidylserine decarboxylase proenzyme, mitochondrial (EC 4.1.1.65) [Cleaved into: Phosphatidylserine decarboxylase beta chain; Phosphatidylserine decarboxylase alpha chain], PISD

Immunogen:

Recombinant Human Phosphatidylserine decarboxylase proenzyme protein (1-375AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A