

PACO41574

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, nucleus, protein tag, cellular protein modification process, mRNA splicing, via spliceosome, positive regulation of protein targeting to mitochondrion.

**Gene ID:**

UBL5

**Uniprot**

Q9BZL1

**Synonyms:**

Ubiquitin-like protein 5, UBL5

**Immunogen:**

Recombinant Human Ubiquitin-like protein 5 protein (1-73AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A