

PACO41586

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

apical plasma membrane, extracellular exosome, antigen binding, antigen transcytosis by M cells in mucosal-associated lymphoid tissue.

Gene ID:

Gp2

Uniprot

Q9D733

Synonyms:

Pancreatic secretory granule membrane major glycoprotein GP2 (Pancreatic zymogen granule membrane protein GP-2), Gp2

Immunogen:

Recombinant Mouse Pancreatic secretory granule membrane major glycoprotein GP2 protein (22-531AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A