

# NIF3L1 Antibody, HRP conjugated



PACO41591

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May function as a transcriptional corepressor through its interaction with COPS2, negatively regulating the expression of genes involved in neuronal differentiation.

**Gene ID:**

NIF3L1

**Uniprot**

Q9GZT8

**Synonyms:**

NIF3-like protein 1 (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 1 protein), NIF3L1, ALS2CR1

**Immunogen:**

Recombinant Human NIF3-like protein 1 protein (1-377AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A