

PACO41606

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

CCDC68

Uniprot

Q9H2F9

Synonyms:

Coiled-coil domain-containing protein 68 (Cutaneous T-cell lymphoma-associated antigen se57-1) (CTCL-associated antigen se57-1), CCDC68


Immunogen:

Recombinant Human Coiled-coil domain-containing protein 68 protein (131-335AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO41606 at dilution of 1:100.