

C17orf75 Antibody, Biotin conjugated



PACO41665

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May have a role in spermatogenesis.

Gene ID:

C17orf75

Uniprot

Q9HAS0

Synonyms:

Protein Njmu-R1, C17orf75

Immunogen:

Recombinant Human Protein Njmu-R1 protein (201-396AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A