

# LINC00846 Antibody



PACO41702

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LINC00846

**Uniprot**

Q9NV44

**Synonyms:**

Putative uncharacterized protein encoded by LINC00846, LINC00846, C21orf77

**Immunogen:**

Recombinant Human Putative uncharacterized protein encoded by LINC00846 protein (1-126AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A