

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Mitochondrial polymerase processivity subunit. Stimulates the polymerase and exonuclease activities, and increases the processivity of the enzyme. Binds to ss-DNA.

**Gene ID:**

POLG2

**Uniprot**

Q9UHN1

**Synonyms:**

DNA polymerase subunit gamma-2, mitochondrial (EC 2.7.7.7) (DNA polymerase gamma accessory 55 kDa subunit) (p55) (Mitochondrial DNA polymerase accessory subunit) (MtPolB) (PolG-beta), POLG2, MTPOLB

**Immunogen:**

Recombinant Human DNA polymerase subunit gamma-2, mitochondrial protein (1-485AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A