

POLG2 Antibody, Biotin conjugated



PACO41785

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Mitochondrial polymerase processivity subunit. Stimulates the polymerase and exonuclease activities, and increases the processivity of the enzyme. Binds to ss-DNA.

Gene ID:

POLG2

Uniprot

Q9UHN1

Synonyms:

DNA polymerase subunit gamma-2, mitochondrial (EC 2.7.7.7) (DNA polymerase gamma accessory 55 kDa subunit) (p55) (Mitochondrial DNA polymerase accessory subunit) (MtPolB) (PolG-beta), POLG2, MTPOLB

Immunogen:

Recombinant Human DNA polymerase subunit gamma-2, mitochondrial protein (1-485AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A