

SERTAD1 Antibody, HRP conjugated



PACO41795

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts at E2F-responsive promoters as coregulator to integrate signals provided by PHD- and/or bromodomain-containing transcription factors. Stimulates E2F1/TFDP1 transcriptional activity. Renders the activity of cyclin D1/CDK4 resistant to the inhibitory effects of CDKN2A/p16INK4A.

Gene ID:

SERTAD1

Uniprot

Q9UHV2

Synonyms:

SERTA domain-containing protein 1 (CDK4-binding protein p34SEI1) (SEI-1) (Transcriptional regulator interacting with the PHD-bromodomain 1) (TRIP-Br1), SERTAD1, SEI1

Immunogen:

Recombinant Human SERTA domain-containing protein 1 protein (1-236AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A