

GSDMC Antibody, HRP conjugated



PACO41863

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytoplasm, microtubule organizing center, mitochondrion.

Gene ID:

GSDMC

Uniprot

Q9BYG8

Synonyms:

Gasdermin-C (Melanoma-derived leucine zipper-containing extranuclear factor), GSDMC, MLZE

Immunogen:

Recombinant Human Gasdermin-C protein (1-508AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A