

CLPB Antibody, Biotin conjugated



PACO41885

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May function as a regulatory ATPase and be related to secretion/protein trafficking process.

Gene ID:

CLPB

Uniprot

Q9H078

Synonyms:

Caseinolytic peptidase B protein homolog (EC 3.6.1.3) (Suppressor of potassium transport defect 3), CLPB, HSP78 SKD3

Immunogen:

Recombinant Human Caseinolytic peptidase B protein homolog protein (1-350AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A