

# KLHL42 Antibody, FITC conjugated



PACO42028

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex required for mitotic progression and cytokinesis. The BCR(KLHL42) E3 ubiquitin ligase complex mediates the ubiquitination and subsequent degradation of KATNA1. Involved in microtubule dynamics throughout mitosis.

**Gene ID:**

KLHL42

**Uniprot**

Q9P2K6

**Synonyms:**

Kelch-like protein 42 (Cullin-3-binding protein 9) (Ctb9) (Kelch domain-containing protein 5), KLHL42, KIAA1340 KLHDC5

**Immunogen:**

Recombinant Human Kelch-like protein 42 protein (1-410AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A