

STOML1 Antibody



PACO42034

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:**Gene ID:**

STOML1

Uniprot

Q9UBI4

Synonyms:

Stomatin-like protein 1 (SLP-1) (EPB72-like protein 1) (Protein unc-24 homolog) (Stomatin-related protein) (STORP), STOML1, SLP1 UNC24


Immunogen:

Recombinant Human Stomatin-like protein 1 protein (79-290AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of A549 cells using PACO42034 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).