

# SLC23A1 Antibody, HRP conjugated



PACO42055

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Sodium/ascorbate cotransporter. Mediates electrogenic uptake of vitamin C, with a stoichiometry of 2 Na(+) for each ascorbate.

**Gene ID:**

SLC23A1

**Uniprot**

Q9UHI7

**Synonyms:**

Solute carrier family 23 member 1 (Na(+)/L-ascorbic acid, transporter 1) (Sodium-dependent vitamin C transporter 1) (hSVCT1) (Yolk sac permease-like molecule 3), SLC23A1, SVCT1 YSPL3

**Immunogen:**

Recombinant Human Solute carrier family 23 member 1 protein (173-259AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A