

SLC23A1 Antibody, Biotin conjugated



PACO42057

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Sodium/ascorbate cotransporter. Mediates electrogenic uptake of vitamin C, with a stoichiometry of 2 Na(+) for each ascorbate.

Gene ID:

SLC23A1

Uniprot

Q9UHI7

Synonyms:

Solute carrier family 23 member 1 (Na(+)/L-ascorbic acid, transporter 1) (Sodium-dependent vitamin C transporter 1) (hSVCT1) (Yolk sac permease-like molecule 3), SLC23A1, SVCT1 YSPL3

Immunogen:

Recombinant Human Solute carrier family 23 member 1 protein (173-259AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A