

C11orf84 Antibody, HRP conjugated



PACO42099

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

C11orf84

Uniprot

Q9BUA3

Synonyms:

SPIN1-docking protein (SPIN-DOC), C11orf84, SPINDOC

Immunogen:

Recombinant Human Uncharacterized protein C11orf84 protein (82-381AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A