

TRIM15 Antibody



PACO42142

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

innate immune response, mesodermal cell fate determination, negative regulation of intracellular transport of viral material, negative regulation of viral release from host cell, positive regulation of NF-kappaB transcription factor activity, positive regulation of RIG-I signaling pathway, positive regulation of sequence-specific DNA binding transcription factor activity, positive regulation of type I interferon production.

Gene ID:

TRIM15

Uniprot

Q9C019

Synonyms:

Tripartite motif-containing protein 15 (RING finger protein 93) (Zinc finger protein 178) (Zinc finger protein B7), TRIM15, RNF93 ZNF178 ZNFB7


Immunogen:

Recombinant Human Tripartite motif-containing protein 15 protein (122-387AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TRIM15 antibody at 2 μ g/ml. Lane 1: Jurkat whole cell lysate. Lane 2: HeLa whole cell lysate. Lane 3: HepG2 whole cell lysate. Lane 4: A549 whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 53, 13 kDa. Observed band size: 53 kDa.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO42142 at dilution of 1:100.