

ROPN1 Antibody



PACO42198

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytoplasm, nucleus, small GTPase mediated signal transduction.

Gene ID:

ROPN1

Uniprot

Q9HAT0

Synonyms:

Ropporin-1A (Cancer/testis antigen 91) (CT91) (Rhophilin-associated protein 1A), ROPN1, ROPN1A


Immunogen:

Recombinant Human Ropporin-1A protein (1-212AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO42198 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO42198 at dilution of 1:100.