

# LRRC4C Antibody, FITC conjugated



PACO42212

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May promote neurite outgrowth of developing thalamic neurons.

**Gene ID:**

LRRC4C

**Uniprot**

Q9HCJ2

**Synonyms:**

Leucine-rich repeat-containing protein 4C (Netrin-G1 ligand) (NGL-1), LRRC4C, KIAA1580 NGL1

**Immunogen:**

Recombinant Human Leucine-rich repeat-containing protein 4C protein (201-500AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A