

LRRC4C Antibody, Biotin conjugated



PACO42213

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May promote neurite outgrowth of developing thalamic neurons.

Gene ID:

LRRC4C

Uniprot

Q9HCJ2

Synonyms:

Leucine-rich repeat-containing protein 4C (Netrin-G1 ligand) (NGL-1), LRRC4C, KIAA1580 NGL1

Immunogen:

Recombinant Human Leucine-rich repeat-containing protein 4C protein (201-500AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A