

ABCF3 Antibody



PACO42246

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Displays an antiviral effect against flaviviruses such as west Nile virus (WNV) in the presence of OAS1B.

Gene ID:

ABCF3

Uniprot

Q9NUQ8

Synonyms:

ATP-binding cassette sub-family F member 3, ABCF3


Immunogen:

Recombinant Human ATP-binding cassette sub-family F member 3 protein (201-500AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO42246 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO42246 at dilution of 1:100.