

# FBXO28 Antibody, HRP conjugated



PACO42259

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probably recognizes and binds to some phosphorylated proteins and promotes their ubiquitination and degradation.

**Gene ID:**

FBXO28

**Uniprot**

Q9NVF7

**Synonyms:**

F-box only protein 28, FBXO28, CENP-30 KIAA0483

**Immunogen:**

Recombinant Human F-box only protein 28 protein (251-368AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A