

FAIM Antibody, Biotin conjugated



PACO42265

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role as an inducible effector molecule that mediates Fas resistance produced by surface Ig engagement in B cells.

Gene ID:

FAIM

Uniprot

Q9NVQ4

Synonyms:

Fas apoptotic inhibitory molecule 1, FAIM, FAIM1

Immunogen:

Recombinant Human Fas apoptotic inhibitory molecule 1 protein (1-179AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A