

# PAQR5 Antibody, HRP conjugated



PACO42283

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Steroid membrane receptor. Binds progesterone. May be involved in oocyte maturation.

**Gene ID:**

PAQR5

**Uniprot**

Q9NXX6

**Synonyms:**

Membrane progesterin receptor gamma (mPR gamma) (Membrane progesterone P4 receptor gamma) (Membrane progesterone receptor gamma) (Progesterone and adipoQ receptor family member 5) (Progesterin and adipoQ receptor family member 5) (Progesterin and adipoQ receptor family member V), PAQR5, MPRG

**Immunogen:**

Recombinant Human Membrane progesterin receptor &gamma; protein (1-51AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A