

PAQR5 Antibody, Biotin conjugated



PACO42285

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Steroid membrane receptor. Binds progesterone. May be involved in oocyte maturation.

Gene ID:

PAQR5

Uniprot

Q9NXX6

Synonyms:

Membrane progesterin receptor gamma (mPR gamma) (Membrane progesterone P4 receptor gamma) (Membrane progesterone receptor gamma) (Progesterone and adipoQ receptor family member 5) (Progesterone and adipoQ receptor family member 5) (Progesterone and adipoQ receptor family member V), PAQR5, MPRG

Immunogen:

Recombinant Human Membrane progesterin receptor γ protein (1-51AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A