

# DKK3 Antibody, Biotin conjugated



PACO42321

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Antagonizes canonical Wnt signaling by inhibiting LRP5/6 interaction with Wnt and by forming a ternary complex with the transmembrane protein KREMEN that promotes internalization of LRP5/6. DKKs play an important role in vertebrate development, where they locally inhibit Wnt regulated processes such as antero-posterior axial patterning, limb development, somitogenesis and eye formation. In the adult, Dkks are implicated in bone formation and bone disease, cancer and Alzheimer disease.

**Gene ID:**

DKK3

**Uniprot**

Q9UBP4

**Synonyms:**

Dickkopf-related protein 3 (Dickkopf-3) (Dkk-3) (hDkk-3), DKK3, REIC

**Immunogen:**

Recombinant Human Dickkopf-related protein 3 protein (22-350AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A