

PACO42322

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Binds to nonmethylated 5'-d(CGG)(n)-3' trinucleotide repeats in the FMR1 promoter. May play a role in regulating FMR1 promoter.

Gene ID:

CGGBP1

Uniprot

Q9UFW8

Synonyms:

CGG triplet repeat-binding protein 1 (CGG-binding protein 1) (20 kDa CGG-binding protein) (p20-CGGBP DNA-binding protein), CGGBP1, CGGBP


Immunogen:

Recombinant Human CGG triplet repeat-binding protein 1 protein (1-167AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO42322 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO42322 at dilution of 1:100.