

APPL1 Antibody, HRP conjugated



PACO42351

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for the regulation of cell proliferation in response to extracellular signals from an early endosomal compartment. Links Rab5 to nuclear signal transduction.

Gene ID:

APPL1

Uniprot

Q9UKG1

Synonyms:

DCC-interacting protein 13-alpha (Dip13-alpha) (Adapter protein containing PH domain, PTB domain and leucine zipper motif 1), APPL1, APPL DIP13A KIAA1428

Immunogen:

Recombinant Human DCC-interacting protein 13-α protein (214-495AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A