

TRAK1 Antibody, Biotin conjugated



PACO42369

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the regulation of endosome-to-lysosome trafficking, including endocytic trafficking of EGF-EGFR complexes and GABA-A receptors.

Gene ID:

TRAK1

Uniprot

Q9UPV9

Synonyms:

Trafficking kinesin-binding protein 1 (106 kDa O-GlcNAc transferase-interacting protein), TRAK1, KIAA1042 OIP106

Immunogen:

Recombinant Human Trafficking kinesin-binding protein 1 protein (101-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A