

HOOK1 Antibody, HRP conjugated



PACO42371

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for spermatid differentiation. Probably involved in the positioning of the microtubules of the manchette and the flagellum in relation to the membrane skeleton. Component of the FTS/Hook/FHIP complex (FHF complex). The FHF complex may function to promote vesicle trafficking and/or fusion via the homotypic vesicular protein sorting complex (the HOPS complex).

Gene ID:

HOOK1

Uniprot

Q9UJC3

Synonyms:

Protein Hook homolog 1 (h-hook1) (hHK1), HOOK1

Immunogen:

Recombinant Human Protein Hook homolog 1 protein (301-600AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A