

ANGPTL2 Antibody



PACO42378

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Induces sprouting in endothelial cells through an autocrine and paracrine action.

Gene ID:

ANGPTL2

Uniprot

Q9UKU9

Synonyms:

Angiotensin-related protein 2 (Angiotensin-like protein 2), ANGPTL2, ARP2


Immunogen:

Recombinant Human Angiotensin-related protein 2 protein (23-493AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO42378 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO42378 at dilution of 1:100.